
RC-CORNET® TRACHEOSET

For using the RC-Cornet® with tracheostomized patients

Each RC-Cornet® Tracheoset includes:

- RC-One way valve
- RC-Cornet® adapter with mouthpiece
- RC-Cornet® adapter to trachea cannula

Please read the RC-Cornet® Tracheoset instructions for use carefully before first use. Please also follow the instructions of the RC-Cornet®.

CEGLA
MEDICAL TECHNOLOGY





Manufacturer:

R. Cegla GmbH & Co. KG
Horresser Berg 1
56410 Montabaur
Germany
Tel +49 2602 9213-0
info@cegla.de
www.cegla.de

- not autoclavable
- not boilable
- not microwaveable

QUICK START GUIDE

PURPOSE FOR USE

The RC-Cornet® Tracheoset allows tracheostomy patients who can breathe spontaneously (even if only for a short period) respiratory physiotherapy with the RC-Cornet®. The RC-Cornet® frees the respiratory tracts effectively from thick bronchial mucus, reduces shortness of breath and coughing.

INSTRUCTIONS FOR USE

I Connect the RC-Cornet® Tracheoset with the RC-Cornet® by sliding the RC-Cornet® onto the adapter support. Push firmly ensuring a tight connection.

II Choose the mouthpiece setting for the desired pressure profile. Pressure is increased as you adjust the settings (0 = Lower, 4 = Higher).

III Connect the RC-Cornet® Tracheoset firmly to the cannula (tubes). When starting the therapy the mouthpiece may be pulled 2 cm out of the curved tube.

IV After a steady breathing phase with a subsequent short break, breathe through the RC-Cornet® Tracheoset in and out. While exhaling, a positive pressure is built in the bronchial system and vibrations are generated. After several breaths, rest as required then repeat. Once you achieved more secretion, please remove the system of the cannula.

V The recommended usage is 2–3 times daily for approx. 2–5 minutes. Patients affected by severe mucus production may require more frequent use.

CLEANING INSTRUCTIONS

The RC-Cornet® Tracheoset is designed for single patient use. It is not autoclavable. Please only clean in detergent. Please substitute the device after 8 days.

The RC-Cornet® can be clinically sterilized in an autoclave for multiple patient use. The instruction for multiple patient use can be found at www.cegla.de.

Benefits of the RC-Cornet® Tracheoset:

- loosens thick bronchial mucus
- widens the respiratory tracts
- improves secretion management
- supports the weaning process (weaning from the respiratory protective device)
- shortens weaning time
- strengthens respiratory muscles
- reduces atelectasis

NOTES

Please also follow the contraindications of the RC-Cornet®.

The RC-Cornet® Tracheoset should only be used with clinical supervision.

The RC-Cornet® Tracheoset can be directly attached to the tracheal cannula.

The one way valve must not be detached from the adapter due to the danger of swallowing. In the event the one way valve becomes detached, then the RC-Cornet® Tracheoset should be replaced.

article-no. 270EN

Replacement parts are available in pharmacies, medical supply stores or at www.cegla-shop.com

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