



RC CLEAN

Professional Clean, with Steam – simple, fast and safe

The RC-Clean cleaning bag allows for an efficient and simple hygienic processing of medical products including baby items in the microwave.

CEGLA
MEDICAL TECHNOLOGY



Cleaning bag for the microwave

Whether in the hospital, at a GP practice, at home or when travelling – the RC-Clean is simple, fast, safe and effective. Microbiological examinations confirm the hygienic effectiveness of the RC-Clean cleaning bag (test report available upon request). All medical products which can be sterilized in boiling water and are microwave safe, such as the RC-Cornet® and the RC-Cornet® N, are suitable for use with the RC-Clean. No metal can be used with the product and it is not working in an oven.

Product Details

The package includes 5 RC-Clean cleaning bags.
The cleaning bag can be reused up to 20 times. Therefore each package offers up to 100 cleaning cycles.

Article-no. 500EN | PPN 111 075 103 265



Available at pharmacies or at:
www.cegla-shop.com

CEGLA
MEDICAL TECHNOLOGY

R. Cegla GmbH & Co. KG
Horresser Berg 1
56410 Montabaur
Germany
Tel +49 2602 92 13-0
info@cegla.de
www.cegla.de